

XTCav processing on psana

The script used to create the data files is scripts/xtcav_process.py

To execute connect to pslogin, then ssh psana, then e.g.

```
bsub -J XCAV.186 -q psanaq -o XCAV.186.log -n 10 mpirun python
xtcav_process.py 186
```

or you can launch multiple with this one-liner:

```
for r in {1..186}; do bsub -J XT.$r -q psanaq -o
/path/to/logfolder/XCAV.$r.log -n 2 mpirun python xtcav_process.py $r;
done
```

It is here for preservation:

```
import sys,os

from scipy.signal import argrelemax
from sklearn.cluster import KMeans
import numpy as np

import psana
from xtcav2.LasingOnCharacterization import LasingOnCharacterization
xx = LasingOnCharacterization()

#####
run = int( sys.argv[1] )
exp = "cxilu5617"
outdir = "/reg/d/psdm/cxi/cxilu5617/scratch/xtcav/batch.data.with.trace"
power_max = 10000 # sometimes power analysis diverges, suppress this
case
max_order = 3 # parameter describing local maxima neighborhood in power
spectrum, increase to lower the local max detected
#####

if not os.path.exists(outdir):
    os.makedirs( outdir)

outname = os.path.join( outdir, "run%d_xtcav.h5"%run)

def smooth(x, beta=10.0, window_size=11):
    """a nice recipe"""
    if window_size % 2 == 0:
        window_size += 1
```

```

s = np.r_[x[window_size-1:0:-1], x, x[-1:-window_size:-1]]
w = np.kaiser(window_size, beta)
y = np.convolve( w/w.sum(), s, mode='valid' )
b = (window_size-1) / 2
smoothed = y[b:len(y)-b]

return smoothed

ds = psana.MPIDataSource("exp=%s:run=%d:smd"%(exp,run))
smldata = ds.small_data(outname) # , gather_interval=100)

# for finding the max position(s) in power spectrum
# we use Kmeans
k_for_max = KMeans(n_clusters=2)

events = ds.events()
####
# first iterate and figure out how long the traces are
for i_ev, ev in enumerate(events):
    xx.processEvent(ev)
    pulse_t, pulse_power = xx.xRayPower()
    if pulse_power is not None:
        Npts = pulse_power[0].shape[0]
        break
    else:
        print "skipping event %d"%i_ev

#events = ds.events()
dummie = np.zeros( Npts)

for i_ev, ev in enumerate( events):
    if ev is None:
        continue
    xx.processEvent(ev)

    pulse_t, pulse_power = xx.xRayPower()
    if pulse_t is None or pulse_power is None:
        print "Power is None %d"%i_ev
        first_pulse_power = -1 #np.nan
        second_pulse_power = -1 #np.nan
        pulse_separation = -1 #np.nan
        first_pulse_time = -1 #np.nan
        second_pulse_time = -1 #np.nan
        first_pulse_max_position = -1 #np.nan
        second_pulse_max_position = -1 #np.nan
        is_separated = -1 #np.nan
        pulse_t = dummie
        pulse_power = dummie
        smooth_power = dummie
    else:

```

```

pulse_t = pulse_t[0]

# filter/cleanup the power trace
pulse_power = np.nan_to_num( pulse_power[0] )
pulse_power[ pulse_power > power_max] = 0

# work with a smoothed power trace
smooth_power = smooth( pulse_power)

# finds all local max in power trace
maxpos = argrelmax( smooth_power, order=max_order)[0]

# handle use cases
if len( maxpos) ==0 :
    first_pulse_power = -1# np.nan
    second_pulse_power = -1# np.nan
    pulse_separation = -1#np.nan
    first_pulse_time = -1#np.nan
    second_pulse_time = -1#np.nan
    first_pulse_max_position = -1#np.nan
    second_pulse_max_position = -1#np.nan
    is_separated = -1#np.nan
    print "No max found baby %d"%i_ev
elif len( maxpos) == 1:
    first_pulse_power = smooth_power[maxpos[0]]
    second_pulse_power = -1# np.nan
    first_pulse_time = pulse_t[ maxpos[0] ]
    second_pulse_time = -1#np.nan
    pulse_separation = -1#np.nan

    first_pulse_max_position = maxpos[0]
    second_pulse_max_position = -1#np.nan
    is_separated = 0 #False
    print "I see one peak in the power spectrum babe; %d"%i_ev
else:
    k_for_max.fit( maxpos[:,None] )
    maxpos1 = maxpos[k_for_max.labels_==0]
    maxpos2 = maxpos[k_for_max.labels_==1]

    maxVals = smooth_power[maxpos]
    maxVals1 = maxVals[k_for_max.labels_==0]
    maxVals2 = maxVals[k_for_max.labels_==1]

# the true peak is the max position within each cluster
X = maxpos1[ np.argmax( maxVals1 )]
X2 = maxpos2[ np.argmax( maxVals2 )]
X,X2 = np.sort( [X,X2])

first_pulse_power = smooth_power[X]
second_pulse_power = smooth_power[X2]

```

```
first_pulse_time = pulse_t[ X ]
second_pulse_time = pulse_t[X2]
pulse_separation = pulse_t[X2] - pulse_t[X]
first_pulse_max_position = X
second_pulse_max_position = X2
is_separated = 1 #True
print "Looks like 2 peaks to me.. %d"%i_ev
smldata.event( first_pulse_power = first_pulse_power,
second_pulse_power = second_pulse_power,
first_pulse_time = first_pulse_time,
second_pulse_time = second_pulse_time,
pulse_separation=pulse_separation,
first_pulse_maxpos = first_pulse_max_position,
second_pulse_maxpos = second_pulse_max_position,
is_separated=is_separated,
power_trace = pulse_power,
smoothed_power_trace = smooth_power, pulse_t=pulse_t)
```

```
#print "processed %d events"%i_ev
smldata.save()
```

The h5 file produced has a lot of useful info, here is a brief example of how to access the data inside:

```
import h5py
from pylab import *

f =
h5py.File('/reg/d/psdm/cxi/cxilu5617/scratch/xtcav/batch.data.with.trace
/run186_xtcav.h5','r')

# display some traces
plot(f["pulse_t"][0], f["power_trace"][0] , 'o', color='C0')
plot(f["pulse_t"][0], f["smoothed_power_trace"][0] , '--', color='C0')
plot(f["pulse_t"][100], f["power_trace"][100] , 'o', color='C1')
plot(f["pulse_t"][100], f["smoothed_power_trace"][100] , '--',
color='C1')
legend( ("shot 1", "smoothed shot 1", "shot 100 ", "smoothed shot 100"))
xlabel("time (fs)")
ylabel("power")

# get information on the shots..
info0 = 0, f['first_pulse_power'][0], f["second_pulse_power"][0],
f["pulse_separation"][0]
info100 =100, f['first_pulse_power'][100], f["second_pulse_power"][100],
f["pulse_separation"][100]
info_s = "event %d, power1: %.2f GW, power2: %.2f GW, separation: %.2f
fs"

title_s = "xtcav traces LU56 from run 186\n" +info_s%info0 + "\n" +
info_s%info100
title(title_s)

show()
```

If you run the above code on psana it should produce the following plot

xtcav traces LU56 from run 186

event 0, power1: 31.36 GW, power2: 25.08 GW, separation: 58.38 fs

event 100, power1: 21.26 GW, power2: 31.25 GW, separation: 46.99 fs

